OPERATION MANUAL

Rx Only

MADE IN JAPAN

OM-DT0919EN 001

Stainless Turbine

S-Max M

S-Max pico

Pana-Max2

Pana-Max PLUS2

1 User and Indications for Use

User: Qualified professionals (dentist)

Indications for Use:

The Stainless Turbine is intended for the following applications: Caries removal, Cavity and crown preparation, Removal of dental restorations (fillings and prostheses), Finishing of teeth and dental restorations (preparation/adjustment).

Patient population: Persons over 2 years old with teeth or artificial teeth.

*** See the “12–2 Symbol” section on page 53 for reference.**

2 Precautions for Handling and Operation

- Please read these precautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
▲ WARNING	Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.

⚠ CAUTION	Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.
NOTICE	General product specification information highlighted to avoid product malfunction and performance reduction.

⚠ WARNING

- Depressing the push button while the handpiece is in operation may cause overheating, serious technical damage and possible premature handpiece failure. During operation avoid contact with any oral tissue that may cause the push button to be depressed while the handpiece is in operation.

⚠ WARNING

- The handpiece is delivered in a non-sterile condition and must be cleaned and sterilized prior to its first use and after each patient use.
- Clean and lubricate the handpiece immediately (within 1 hour) after each treatment to remove residue.
Failure to properly maintain the handpiece may cause infection, product failure or overheating leading to burn injuries. (Refer to “5 Post-use Maintenance”)
- If blood infiltrates inside a handpiece, clean and lubricate the handpiece immediately using the PANA SPRAY Plus or clean the handpiece using a washer-disinfector.

WARNING

- Keep any debris or other foreign materials away from inside the cartridge or the handpiece. Foreign materials remaining inside may lead to overheating, causing burn injuries or other accidents.
- Ball bearings are wear items and may require replacement. Be sure to inspect the handpiece before use (Refer to “4–4 Checking the Handpiece Before Each Use”). If abnormal vibration or noise are found, the bearings or other internal parts may be worn out or damaged. Damaged parts may lead to handpiece overheating causing burn injuries. If overheating is suspected, stop using handpiece immediately and contact your Authorized NSK Dealer.

CAUTION

- Read this operation manual before use to fully understand the product functions and file for future reference.
- When operating the product always consider the safety of the patient.
- Users are responsible for the operational control, maintenance and continual inspection of this product.
- Prior to clinical use, inspect the handpiece. Check for vibration, noise and overheating. If any abnormalities are found, stop using the handpiece immediately and contact your Authorized NSK Dealer. (Refer to “4–4 Checking the Handpiece Before Each Use”)
- Do not disassemble or alter the handpiece except as recommended by NSK in this operation manual.

CAUTION

- Do not allow any impact on to the product. Do not drop the product. Deformation may cause the handpiece to fail during use.
- Operators and all others in the area must wear eye protection and a mask when operating this handpiece.
- An accidental symptom, such as subcutaneous emphysema may be developed in the use of this product. Do not direct air exhausted from the handpiece to the bleeding area (wound) or into periodontal pockets.
- Should the handpiece function abnormally during use, stop using the handpiece immediately and contact your Authorized NSK Dealer.

CAUTION

- The handpiece is delivered in a non-sterile condition and must be cleaned and sterilized prior to use.
- Do not use the following fluids to wipe, immerse or clean the product; strong/super acid water, strong acid/alkaline chemicals, solutions containing chlorine, solvents such as benzine or thinner.
- Perform regular function and maintenance checks. (Refer to “10–4 Periodical Maintenance Checks”)
- If the handpiece has not been used for a long period, rotate the handpiece and check for noise, vibration and overheating before use.
- To avoid clinical downtime it is recommended that a spare be kept on hand in case of a breakdown during treatment.

⚠ CAUTION

- Do not adverse effect the pulp and use it intermittently with a feather touch under water injection.
- If water leaks from the tip of the handpiece even after the handpiece operation is stopped, replace the Anti-Retraction Valve. (Refer to “10–3 Replacing the Anti-Retraction Valve” (S-Max pico BLED only)) For products other than S-Max pico BLED, replace the Anti-Retraction Valve according to the coupling’s instruction manual.

M800WL / M900WL / S-Max pico BLED / S-Max pico WLED

- Avoid continual eye contact with the LED light.

⚠ CAUTION

- This product is classified as Class 1 LED Product.
- Use a power source which meets the following requirements.
 - 1 The electricity supply of the power source is below 15W both under normal and single-failure conditions.
 - 2 The power source uses a SELV circuit for electricity supply.
 - 3 The output voltage of the power source is within the range recommended by the manufacturer of this product.
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

NOTICE

- Please use the handpiece in a place with illuminance of 500lx or more.
- Patients may feel pain due to vibration and heat transferred to the pulp during treatment.
- Repeated use may wear the inside of the head cap, which may make it difficult to push the push button. In this case, please replace the head cap.

3 Package Contents

No.	Part Name	Quantity	Remarks
1	Spray Nozzle	1	Excluding Pana-Max2 M4 / PAP2-MU M4 / PAP2-SU M4
2	Cleaning Wire	1	-
3	PTL O-ring Set	1	Included in the package for M800L / M800 / M900L / M900 / S-Max pico
4	Head Cap Wrench	1	-

4 Preparation for Use

4-1 Setting of Air & Water Supply Pressure

Measure the supply pressure at the handpiece/hose connection point and set the pressure to the value specified in “12-1 Specifications”. (Fig. 1).

For NSK multi gauge information refer to Optional Spare Parts List.


Fig. 1

Measurement may not be possible at the handpiece/hose connection point depending on the specifications of the hose for the dental unit. Contact the manufacturer of the dental unit for details.

⚠ WARNING

- Do not exceed the optimum pressure specified in “12-1 Specifications”.

⚠ CAUTION

- Do not use air contaminated by dust, moisture or oil.

4-2 Connecting and Disconnecting the Handpiece

4-2-1 Connecting

M800L / M800 / M900L / M900 / S-Max pico

- 1 Insert the handpiece into the coupling (Fig. 2).
- 2 Make sure the handpiece is fully seated in the coupling (Fig. 3).


Fig. 2


Fig. 3

Pana-Max2 M4 / PAP2-MU M4 / PAP2-SU M4

- 1 Insert the handpiece correctly into the hose connector and tighten the hose nut (Fig. 4).
- 2 Make sure the handpiece is connected firmly to the hose (Fig. 5).


Fig. 4


Fig. 5

4-2-2 Disconnecting

M800L / M800 / M900L / M900 / S-Max pico

Pull back the retention lock ring and remove the handpiece from the coupling (Fig. 6).

Pana-Max2 M4 / PAP2-MU M4 / PAP2-SU M4

Loosen the hose nut and remove the handpiece from the hose (Fig. 7).


Fig. 6


Fig. 7

4-2-3 Connecting and Disconnecting The other handpieces

M800KL / M800WL / M900KL / M900K / M900WL / M900W / S-Max pico KL / S-Max pico SL / S-Max pico BLED / S-Max pico WLED

For the handpieces listed above, connecting and disconnecting the handpiece according to the operation manual for coupling to be used.

⚠ CAUTION

- Do not operate the retention lock ring while pressurized. This may cause the handpiece to shoot off of the coupling unexpectedly.

4-3 Mounting and Removing the Bur

4-3-1 Mounting the Bur

- 1 Insert the bur into the chuck (Fig. 8).
- 2 Depress the push button to open the chuck (1) (Fig. 9).
- 3 Insert the bur fully into the chuck until it stops (2) then release the push button (Fig. 9).
- 4 Make sure that the bur is not removed by gently pulling and pushing the bur without depressing the push button (Fig. 10).

NOTICE

- Grip the handpiece while placing your thumb tip on the push button which makes it easier to depress the button.

4-3-2 Removing the Bur

Depress the push button to open the chuck (1) and remove the bur (2) (Fig. 11).


Fig. 8


Fig. 9


Fig. 10


Fig. 11

WARNING

- Always use a clean, sterile bur. Be sure to wear sterile gloves before mounting the bur.
- When changing burs during treatment, remove the bur and wipe around the bur insertion hole and gloves with a dry lint free cloth until the insertion hole is visually clean, then mount a clean, sterile bur. Residual blood or debris may enter the handpiece while changing burs, which could lead to higher risk of cross-contamination.

CAUTION

- Do not use burs with problems listed below as the bur may break, seize up or disengage from the chuck.
 - Bent, deformed, worn, rusted, broken, deficient bur.
 - Bur which is cracked on the edge or axis.
 - Non-ISO standard, or tampered bur.
- Do not use burs other than specified in “12–1 Specifications” as such burs may be accidentally released or break during rotation.
- Insert a sterilized bur and always keep the bur shank clean. When changing burs, make sure that there is no dirt or debris in and around the bur insertion hole of the handpiece. Debris in the chuck could cause infection, heating, poor bur concentricity or low chuck retention force.

CAUTION

- Always insert the bur all the way into the chuck. If insertion is insufficient, premature failure of the bearings or accidental release of the bur may occur.
- Do not exceed Max. Bur Length shown in “12-1 Specifications”.
- Always follow the instructions provided by the bur manufacturer.
- Do not exceed the bur speed recommended by the bur manufacturer.
- Do not mount or remove the bur until the drive air has completely stopped.

CAUTION

- Do not apply excess pressure to the bur as it may break, bend, or become difficult to remove, or the chuck may wear which causes the bur to detach. If the bur cannot be removed, stop using the handpiece and contact your Authorized NSK Dealer.
- Entry of debris into the chuck could cause bur rotation slip and also prevents the bur from being securely located in the chuck.

4-4 Checking the Handpiece Before Each Use

Follow the check below before use. If any abnormalities are found, stop using the handpiece immediately and contact your Authorized NSK Dealer.

- 1 Check the head cap is firmly tightened (Fig. 12).
If it is loose, tighten it with the head cap wrench.
- 2 Mount the bur. (Refer to “4-3 Mounting and Removing the Bur”)
- 3 Rotate the handpiece for about 20 seconds. During rotation, check for abnormalities such as abnormal rotation, vibration, noise. Check the coolant water is flowing properly.
When using an optic handpiece or coupling, make sure that the light is illuminated properly (Fig. 13).

- 4 After the handpiece rotation has completely stopped, touch the handpiece head to confirm the head is NOT heating abnormally (Fig. 14).


Fig. 12

Fig. 13


Fig. 14

CAUTION

- To avoid injury, keep your hands away from the bur during rotation.

5 Post-use Maintenance

After each patient maintain the product as follows

WARNING

- To ensure cleaning and sterilization efficacy, use only the following procedures for reprocessing.
- Follow any additional local directives, standards, and guidelines for cleaning and sterilization.
- The handpiece is delivered in a non-sterile condition and must be cleaned and sterilized prior to its first use and after each patient use.
- Clean and lubricate the handpiece immediately (within 1 hour) after each treatment to remove residue. Failure to properly maintain the handpiece may cause infection, product failure or overheating leading to burn injuries.

5-1 Preparation

⚠ CAUTION

- Do not use the following liquids to wipe, immerse, or clean the product: strong/super acid water, strong acid/alkaline chemicals, solutions containing chlorine, solvents such as benzine or thinner.
- Do not immerse NSK instruments in disinfectant solution or clean in ultrasonic devices.

- 1 Always wear protective gloves, a mask, and protective goggles for safety purposes and to minimize the risk of infection (Fig. 15).
- 2 Remove the bur (Fig. 16).
- 3 Wipe the exterior of the handpiece and coupling clean using a dry cloth* (Fig. 17).

* Lint free cloth is recommended.

Fig. 15


Fig. 16


Fig. 17

- 4 Operate the handpiece at the chair-side for at least 20 seconds to purge fluids from the handpiece (Fig. 18). After the handpiece rotation has completely stopped, touch the handpiece head to confirm the head is NOT heating abnormally.


Fig. 18

⚠ CAUTION

- If the head is heating abnormally, stop using the handpiece immediately and contact your Authorized NSK Dealer.
- 5 Remove the handpiece from the coupling or hose (Fig. 19). Perform maintenance of the coupling according to the operation manual for the coupling.
- 6 Carry the handpiece to the decontamination area.


Fig. 19

5-2 Cleaning and Drying

Perform cleaning and drying using either the manual or automatic method.

Manual Method (Cleaning the Exterior)

Wash the handpiece under running water following the procedures below.

(Water condition: $\leq 38\text{ }^{\circ}\text{C}$, $\geq 3.5\text{ L/min}$; water should be of the same quality as drinking water)

- 1 Clean the external surfaces of the handpiece with a soft bristled toothbrush for 15 seconds or more (Fig. 20).
- 2 Clean around the bur insertion hole for 15 seconds or more, using an interdental brush (wire thickness 0.7mm) that conforms to size 4 in ISO 16409 (Fig. 20).

Under appropriate lighting (natural office lighting, 500 lx or higher), inspect the handpiece for blood or other debris. If any visible debris remains, repeat the process until the handpiece is visually clean.


Fig. 20

Optional

Wipe the exterior of the handpiece clean using a cloth moistened with disinfectant or ethanol with a concentration of 60 - 90 %. (Fig. 21).


Fig. 21

NOTICE

- Wiping is not an alternative method of cleaning. Please do not skip the required cleaning steps which are previously displayed.

When using disinfectant, follow the instructions given by the manufacturer of the disinfectant.

Use a state-sanctioned disinfectant with proven bactericidal, fungicidal, and virucidal properties.

The following disinfectants can be used in the United States and Canada:

CaviCide[®], CaviWipes[®] (manufactured by Metrex)

After cleaning, if moisture is present, wipe it off with dry cloth or blow it off with compressed air (≤ 0.35 MPa) until there is no moisture in the interior and exterior of the handpiece.

When blowing off with compressed air, cover the handpiece with a cloth to prevent the scattering of water. Proceed to "5-3 Cleaning and Lubrication".

Automatic Method (Cleaning and Drying the Exterior and Interior)

This icon denotes that the product can be washed via washer-disinfector

NSK handpieces with this symbol are compatible with medical washer-disinfector.

Due to the variation in cleaning/disinfection cycles and cleaning agents available from the various manufacturers, NSK has validated the automatic processing/reprocessing of the product using Cycle P7-Anesthesia of the Getinge 46 Series washer-disinfector.

Only validated procedures must be used for cleaning and sterilization. When using reprocessing procedures that are different from those described in this manual, those procedures must be validated by the respective practice or hospital using a washer-disinfector that conforms to ISO 15883-1.

NOTICE

- All water is “building supply” (at least drinking water quality) unless otherwise indicated.

- 1 Prepare and pre-cleaning the device as described above in step “5–1 Preparation”.
- 2 Wipe off dirt on a part held by a handpiece holder using a cloth moistened with 60 - 90 % ethanol or disinfectant (CaviWipes®) (Fig. 22).


Fig. 22

When using disinfectant, follow the instructions given by the manufacturer of the disinfectant.

Use a state-sanctioned disinfectant with proven bactericidal, fungicidal, and virucidal properties.

The following disinfectants can be used in the United States and Canada:

CaviCide[®], CaviWipes[®] (manufactured by Metrex)

- 3 Perform main wash using a washer-disinfector that conforms to ISO 15883-1. NSK has validated the following conditions (Cycle P7-Anesthesia of the Getinge 46 Series washer-disinfector).

Pre-Wash 1: Cold water (<110 °F (43 °C)) for 2 min

Pre-Wash 2: None - Not Applicable

Enzyme Wash: Hot tap water and 60ml of heated (95 °F (35 °C)) neutral enzymatic detergent for 3 min at 120 °F (49 °C)

Detergent Wash: None - Not Applicable

Rinse 1: Hot tap water for 2 min

Rinse 2: Hot tap water for 2 min

R0 Final Rinse: Hot tap water for 32 min at 167 °F (75 °C)

Drying: 20 min

Upon removal from the washer-disinfector, if moisture is present, wipe it off with dry cloth or blow it off with compressed air (≤ 0.35 MPa) until there is no moisture in the interior and exterior of the handpiece.

- 4 After cleaning and drying the handpiece, confirm that it is clean under appropriate lighting (500 lx or higher). If dirt still remains, repeat the procedures.
- 5 Proceed to “5–3 Cleaning and Lubrication”.

⚠ CAUTION

- Handpieces must be removed from the washer-disinfector immediately (within 1 hour) after the cleaning and drying cycle is complete to prevent corrosion.
- After washing with washer-disinfector, prior to lubrication, dry the product until all internal moisture is completely removed. Washer-disinfector moisture remaining inside the product could reduce the effect of lubrication and could cause corrosion inside of the product.

5-3 Cleaning and Lubrication

Cleaning and Lubricating the Interior

- ① Attach the spray nozzle included in the nozzle of the PANA SPRAY Plus (Fig. 23).

Only for Pana-Max2 M4 / PAP2-MU M4 / PAP2-SU M4

Attach a Tip Nozzle to the nozzle of the PANA SPRAY Plus (Fig. 24).


Fig. 23


Fig. 24

- 2 Shake the can 3 or 4 times.
- 3 Insert the spray nozzle into the rear of the handpiece while holding the head of the handpiece with a cloth etc.. (Fig. 25).

Only for Pana-Max2 M4 / PAP2-MU M4 / PAP2-SU M4

Insert the Tip Nozzle into the rear of the handpiece while holding the head of the handpiece with a cloth etc.. (Fig. 26, Fig. 27).

- 4 Hold the handpiece, and spray for 2 or 3 seconds until lubricant comes out of the handpiece head.

Repeat lubrication until foreign material stops coming out of the tip.


Fig. 25


Fig. 26

Fig. 27

⚠ CAUTION

- Do not use sprays other than NSK PANA SPRAY Plus. Using sprays other than PANA SPRAY Plus may cause the handpiece to overheat.
- Use only NSK spray nozzle and Tip Nozzle. The effects of other products are not confirmed.
- Hold the spray can upright.
- Firmly hold the handpiece to prevent it slipping when spray pressure is applied.
- Spray lubricant until it expels from the handpiece head.

⚠ CAUTION

- If you wish to purge excessive oil from inside the handpiece, rotate the handpiece for approx 15 seconds without a bur. During the rotation, do not depress the push button, especially when using absorbent cloth to prevent lubricant from scattering. Depressing the push button during rotation may damage the chuck mechanism (Fig. 28).


Fig. 28

NOTICE

- NSK recommends the use of “Spray Mist Absorber” (REF : Y900084) to prevent oil mist expelling out of the handpiece head.

5-4 Packaging, Sterilizing, and Drying

- Use an FDA-approved sterilization pouch that conforms to ISO 11607-1, and seal the pouch. Maximum temperature for NSK instruments is 135 °C (275 °F).

- Perform steam sterilization with the following conditions.

Type	Gravity Displacement	Pre-Vacuum (Dynamic Air Removal)
Temperature	132 °C (270 °F)	132 °C (270 °F)
Full Cycle Time	15 min	4 min
Drying Time	30 min	30 min


Fig. 29

⚠ CAUTION

- Use an FDA-approved steam sterilizer to perform sterilization.
- Follow local rules, regulations, and guidelines regarding the reprocessing of devices.
- Immediately after sterilization, the product is hot and should be handled with caution.
- Place each product in an individual sterilization pouch for steam sterilization. This is to prevent possible discoloration and damage to the product due to chemical residue from other instruments.
- Clean and lubricate the handpiece prior to sterilization. If blood remains on the internal surface it can become clotted and cause product failure.

⚠ CAUTION

- Do not heat or cool the product too quickly. Rapid change in temperature could cause damage to the product.
- Be sure to use sterilizers that can perform sterilization up to 135 °C. In some sterilizers, the chamber temperature may exceed 135 °C. Do not use these sterilizers as failure of the handpiece could occur. Contact the sterilizer manufacturer for detailed information about cycle temperatures.
- Steam sterilization is recommended for the product. The validity of other sterilization methods (such as plasma sterilization or EOG sterilization) is not confirmed.

5-5 Storage

Store the product in a dry, clean location. The product should remain in the sterilization pouch until required for use (Fig. 30).


Fig. 30

⚠ CAUTION

- After the sterilization and drying cycles are complete, remove the handpiece immediately from the sterilizer to store it.
- Store the product in a well ventilated place out of direct sunlight and within the range of temperature, humidity and pressure specified in “12-1 Specifications”.
- Sterilization is not guaranteed after the sterilization retention period specified by the manufacturer and seller of the sterilization pouch has elapsed. If the sterilization retention period has elapsed, perform sterilization again with a new sterilization pouch.

6 Cleaning of Optic Illumination Points

When debris is attached to optic illumination points (Glass rod), wipe clean all the optic illumination points using an alcohol-immersed cotton swab. Remove all debris (Fig. 31).

⚠ CAUTION

- Clean the glass rod after cleaning the handpiece.
- Do not use a sharp tool to clean the glass rod. It could damage the glass and reduce the light transmission.


Fig. 31

7 Cleaning of Spray Ports

- 1 When spray ports are clogged, or spray does not exit evenly from the ports, clean the ports (Fig. 32).


Fig. 32

⚠ CAUTION

- Clean the spray ports while connected to the hose.
- After cleaning the spray ports, perform maintenance from “5-1 Preparation”.
- When cleaning the spray ports during treatment, make sure that adequate cooling water is present outside the patient’s mouth before use.
- Do not forcibly insert the wire into the port. Resultant port damage could cause the spray to be directed away from the bur, and the cooling efficiency reduces.
- Do not blow air into the clean head holes.

- After using the Cleaning Wire, clean the wire for 15 seconds or more by pulling its both ends with your fingers while running water all around the wire (Fig. 33).


Fig. 33

- Under appropriate lighting (500 lx or higher), inspect the Cleaning Wire for blood or debris. If any visible debris remains, repeat the process until they are visually clean.
- If the residual moisture is present after cleaning, wipe it off with dry cloth.
- Sterilize the cleaning wire. (Refer to “5-4 Packaging, Sterilizing, and Drying”.)

8 Cleaning of Spray Nozzle, Tip Nozzle

Clean the spray nozzle and Tip Nozzle under running water according to the following procedures for each use. (Water condition: ≤ 38 °C, ≥ 3.5 L/min; water should be of the same quality as drinking water)

- Clean the entire outer surface of the handpiece with a soft-bristled toothbrush for at least 15 seconds. (Fig. 34). For spray nozzles, clean the knurled part and the O-ring thoroughly. And for Tip Nozzles, clean the triangular part thoroughly.


Fig. 34

- 2 Clean the inside and outside of the spray nozzle and Tip Nozzle with running water.
- 3 Under appropriate lighting (natural office lighting, 500 lx or higher), inspect the spray nozzle and Tip Nozzle for blood or debris. If any visible debris remains, repeat the process until they are visually clean.
- 4 If the residual moisture is present after cleaning, wipe it off with dry cloth.

9 Lubricating the Chuck

Before lubricating the chuck, first clean and dry. (Refer to “5–2 Cleaning and Drying”.)

Manual Method

- 1 Attach the Tip Nozzle to the nozzle of the PANA SPRAY Plus.
- 2 Shake the can 3 or 4 times, then gently press in the push button to directly spray into the bur insertion hole (Fig. 35).


Fig. 35

⚠ CAUTION

- Perform this procedure only for the handpieces that “5–2 Cleaning and Drying” has been completed. If it is performed for the handpieces that “5–2 Cleaning and Drying” is not completed, the inside of the handpiece and the nozzle may be contaminated.
- If the chuck is not regularly Lubricated the chuck grip may be weakened and the bur may be accidentally released during use.

10 Maintenance

10–1 Replacing the Turbine Cartridge

If any of the following symptoms occur, replace the cartridge of the handpiece.

- No rotation, Unstable rotation
- Heat generation
- Bur runout
- Loose bur

- 1 Perform cleaning and sterilizing the handpiece.
- 2 Insert a test bur.
- 3 Locate the correct wrench tool on the head cap then turn the wrench counter clockwise to loosen the cap. Remove the cap (Fig. 36).
- 4 Use the bur to gently lever the entire cartridge out from the head (Fig. 37).


Fig. 36


Fig. 37

- 5 Spray PANA SPRAY Plus into the head for 2 to 3 seconds. Also spray it briefly into the head cap (Fig. 38).
- 6 Wipe the PANA SPRAY Plus lubricant over the head interior.
- 7 Insert the new cartridge into the head by aligning the pin on the cartridge with the slot on the head (Fig. 39). (S-Max pico has no pin and slot.)


Fig. 38


Fig. 39

- ⑧ Attach and tighten the head cap by hand. Then tighten it firmly with the head cap wrench (Fig. 40).
- ⑨ Clean and sterilize the handpiece after replacing the cartridge.
* Refer to Accessories and Spare Parts List to identify the correct parts.


Fig. 40

⚠ CAUTION

- The head cap wrench should only be used on a reprocessed handpiece and outside of the clinical operatory space.
- Use only a genuine NSK cartridge.
- If another cartridge is used NSK cannot guarantee performance and the handpiece warranty would become invalid.
- Always first finger tighten the head cap then secure firmly with the head cap wrench.
- NSK never recommends the disassembly and repair of any NSK cartridge. There is NO EXCEPTION. In such a case a handpiece may perform abnormally (abnormal noise or abnormal vibration). Damage, failure or accidents are outside of our guarantee.

⚠ CAUTION

- Make sure that the O-rings are correctly located on the cartridge (excluding Pana-Max2) (Fig. 41).


Fig. 41

10-2 Replacing the O-rings

(M800L / M800 / M900L / M900 / S-Max pico)

Replace the O-rings if water is present in the exhaust air line. This is an indication of possible water leakage within the coupling. ALWAYS change the complete set of O-rings.

- 1 Loosen and remove the taper ring at the rear of the handpiece (Fig. 42).
 - 2 Gently remove each O-ring by hand (Fig. 43).
 - 3 Insert the complete set of new O-rings in the correct grooves (Fig. 44).
 - 4 Firmly tighten the taper ring (Fig. 45).
- * Refer to Accessories and Spare Parts List to identify the correct parts.


Fig. 42


Fig. 43


Fig. 44


Fig. 45

⚠ CAUTION

- Do not force the new replacement O-ring with excessive force.
- When inserting new O-rings, make sure they are inserted in the correct grooves (Fig. 46).
- Make certain that the taper ring is firmly tightened. If the taper ring is loose water and air leakage could occur.


Fig. 46

10-3 Replacing the Anti-Retraction Valve

If water leaks from the tip of the handpiece even after the handpiece operation is stopped, replace the Anti-Retraction Valve. For products other than the following, replace the Anti-Retraction Valve according to the operation manual for the coupling.

S-Max pico BLED

- 1 Remove the handpiece from the coupling and loosen the taper ring to remove.
- 2 Pull out the joint then pull out the old Anti-Retraction Valve with a thin stick. (Fig. 47).
- 3 Insert new Anti-Retraction Valve securely.
- 4 Mount the joint into the handpiece. Make certain the pipe is aligned. (Fig. 48).
- 5 Tighten the taper ring firmly.


Fig. 47


Fig. 48

⚠ CAUTION

- Make certain that the taper ring is firmly tightened. If the taper ring is loose water and air leakage could occur.

10-4 Periodical Maintenance Checks

Perform periodical maintenance checks every three months, referring to the table below. If any abnormalities are found, contact your Authorized NSK Dealer.

Points to check	Details
Head cap is loose	Check that the head cap is firmly tightened. If it is loose, tighten it with the head cap wrench (Fig. 49).
Rotation	Rotate the handpiece (20 sec.) and check for abnormalities such as abnormal rotation, vibration, noise, and overheating (Fig. 50).
Coolant water	Operate the handpiece and check that the coolant water is flowing through all spray ports (Fig. 51).


Fig. 49


Fig. 50


Fig. 51

11 After-Sales Service

11-1 Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non NSK parts installed. Replacement parts are available for seven years beyond discontinuation of the model. Contact your Authorized NSK Dealer if repairs are necessary.

11-2 Accessories and Spare Parts List

Model	REF	Handpiece Series and Model
Accessories		
Spray Nozzles		
Phatelus Spray Nozzle	P183090	<u>S-Max M</u> M800, M800L, M900, M900L <u>S-Max pico</u> S-Max pico
KV Nozzle	Z090051	<u>S-Max M</u> M800KL, M900K, M900KL <u>S-Max pico</u> S-Max pico KL
SR Nozzle	Z090052	<u>S-Max pico</u> S-Max pico SL
BA Nozzle	Z090053	<u>S-Max pico</u> S-Max pico BLED

Model	REF	Handpiece Series and Model
Accessories		
Spray Nozzles		
WH Nozzle	Z090054	<u>S-Max M</u> M800WL, M900W, M900WL <u>S-Max pico</u> S-Max pico WLED
Cleaning Wire		
Cleaning Wire	Z144101	All models within each handpiece series
Head Cap Wrench		
SX-SU05	20000535	<u>S-Max M</u> M800, M800L, M800KL, M800WL, M900, M900L, M900K, M900KL, M900W, M900WL <u>Pana-Max PLUS2</u> PAP2-MU M4, PAP2-SU M4

Model	REF	Handpiece Series and Model
Accessories		
Head Cap Wrench		
SX-PU05	20002803	<u>S-Max pico</u> S-Max pico, S-Max pico KL, S-Max pico SL, S-Max pico WLED, S-Max pico BLED
PAX2-SU05	P1190765	<u>Pana-Max2</u> Pana-Max2 M4
Spare Parts		
Cartridge		
SX3-MU03	P1280	<u>S-Max M</u> M800, M800L, M800KL, M800WL <u>Pana-Max PLUS2</u> PAP2-MU M4

Model	REF	Handpiece Series and Model
Spare Parts		
Cartridge		
SX3-SU03	P1281	<u>S-Max M</u> M900, M900L, M900K, M900KL, M900W, M900WL <u>Pana-Max PLUS2</u> PAP2-SU M4
SX-PU03	P1140200	<u>S-Max pico</u> S-Max pico, S-Max pico KL, S-Max pico SL, S-Max pico WLED, S-Max pico BLED
PAX2-SU03	P1190	<u>Pana-Max2</u> Pana-Max2 M4

Model	REF	Handpiece Series and Model
Spare Parts		
Head Cap		
SX-MU04	P1001500	<u>S-Max M</u> M800, M800L, M800KL, M800WL <u>Pana-Max PLUS2</u> PAP2-MU M4
SX-SU04 (YOS)	P1027500	<u>S-Max M</u> M900, M900L, M900K, M900KL, M900W, M900WL <u>Pana-Max PLUS2</u> PAP2-SU M4
SX-PU04	P1140500	<u>S-Max pico</u> S-Max pico, S-Max pico KL, S-Max pico SL, S-Max pico WLED, S-Max pico BLED

Model	REF	Handpiece Series and Model
Spare Parts		
Head Cap		
PAX2-SU04	P1190050	<u>Pana-Max2</u> Pana-Max2 M4
O-ring Set		
PTL O-ring Set	Y900580	<u>S-Max M</u> M800, M800L, M900, M900L <u>S-Max pico</u> S-Max pico
Anti-Retraction Valve		
Anti-Retraction Valve	P826042	<u>S-Max pico</u> S-Max pico BLED

11-3 Optional Parts List

Model	REF	Handpiece Series and Model
Multi Gauge		
MG-4H	Z109400	<u>Pana-Max PLUS2</u> PAP2-MU M4, PAP2-SU M4
Nozzle		
Tip Nozzle	Z020010	All models within each handpiece series

11-4 Disposing product

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose of the product for you.

12 Specifications

12-1 Specifications

S-Max M series

Model	M800L	M800	M800KL	M800WL
Coupling Type	NSK PTL		KaVo MULTIflex	W&H Roto Quick
Head Type	Miniature Head			
Rotation Speed	380,000 - 450,000 min ⁻¹			
Bur Type	ISO 1797-1 Type 3 Ø1.59 - 1.60 mm Short Shank Bur / Standard Bur			
Chuck Type	Push Button			
Chucking Length	9 mm			
Max. Bur Length	21 mm			
Max. Working Part Diameter	Ø2 mm			

Model	M800L	M800	M800KL	M800WL
Recommended Drive Air Pressure	0.28 MPa (2.8 kgf/cm ²)			
Operating Drive Air Pressure Range	0.25 - 0.30 MPa (2.5 - 3.0 kgf/cm ²)			
Max. Air Consumption	< 66 NL/min (0.3 ± 0.1 MPa)			
Water Pressure	0.08 - 0.20 MPa (0.8 - 2.0 kgf/cm ²)			
Max. Water Consumption	≥50 mL/min (0.20 MPa)			
Chip Air Pressure	0.10 - 0.20 MPa (1.0 - 2.0 kgf/cm ²)			
Rated Torque	≥0.05 cN·m (Recommended Max Drive Air Pressure)			
Optic	Glass Rod	-	Glass Rod	White LED
Voltage	-			AC/DC 3.3 ± 0.05 V
Active Voltage	-			2.8 - 4.0 V
Consumption Current	-			0.38 A Typical (3.3 V)

Model	M900L	M900	M900KL	M900K	M900WL	M900W
Coupling Type	NSK PTL		KaVo MULTiflex		W&H Roto Quick	
Head Type	Standard Head					
Rotation Speed	325,000 - 430,000 min ⁻¹					
Bur Type	ISO 1797-1 Type 3 Ø1.59 - 1.60 mm Standard Bur					
Chuck Type	Push Button					
Chucking Length	10.5 mm					
Max. Bur Length	25 mm					
Max. Working Part Diameter	Ø2 mm					
Recommended Drive Air Pressure	0.28 MPa (2.8 kgf/cm ²)					
Operating Drive Air Pressure Range	0.25 - 0.30 MPa (2.5 - 3.0 kgf/cm ²)					
Max. Air Consumption	< 66 NL/min (0.3 ± 0.1 MPa)					
Water Pressure	0.08 - 0.20 MPa (0.8 - 2.0 kgf/cm ²)					
Max. Water Consumption	≥50 mL/min (0.20 MPa)					

Model	M900L	M900	M900KL	M900K	M900WL	M900W
Chip Air Pressure	0.10 - 0.20 MPa (1.0 - 2.0 kgf/cm ²)					
Rated Torque	≥0.05 cN·m (Recommended Max Drive Air Pressure)					
Optic	Glass Rod				White LED	-
Voltage	-				AC/DC 3.3 ± 0.05 V	-
Active Voltage	-				2.8 - 4.0 V	-
Consumption Current	-				0.38 A Typical (3.3 V)	-

S-Max pico series

Model	S-Max pico	S-Max pico KL	S-Max pico SL	S-Max pico WLED	S-Max pico BLED
Coupling Type	NSK PTL	KaVo MULTiflex	Sirona Quick Coupling	W&H Roto Quick	Bien-Air Unifix
Head Type	Miniature Head				
Rotation Speed	380,000 - 450,000 min ⁻¹				
Bur Type	ISO 1797-1 Type 3 Ø1.59 - 1.60 mm Short Shank Bur / Super Short Shank Bur				
Chuck Type	Push Button				
Chucking Length	8 mm				
Max. Bur Length	18.5 mm				
Max. Working Part Diameter	Ø2 mm				
Recommended Drive Air Pressure	-				

Model	S-Max pico	S-Max pico KL	S-Max pico SL	S-Max pico WLED	S-Max pico BLED
Operating Drive Air Pressure Range	0.20 - 0.25 MPa (2.0 - 2.5 kgf/cm ²)	0.25 - 0.30 MPa (2.5 - 3.0 kgf/cm ²)			
Max. Air Consumption	45 ± 5 NL/min (0.25 MPa)	45 ± 5 NL/min (0.28 MPa)	45 ± 5 NL/min (0.26 MPa)	45 ± 5 NL/min (0.30 MPa)	45 ± 5 NL/min (0.27 MPa)
Water Pressure	0.05 - 0.20 MPa (0.5 - 2.0 kgf/cm ²)				
Max. Water Consumption	≥50 mL/min (0.20 MPa)				
Chip Air Pressure	0.15 - 0.25 MPa (1.5 - 2.5 kgf/cm ²)				
Rated Torque	≥0.05 cN·m (Recommended Max Drive Air Pressure)				
Optic	Glass Rod			White LED	
Voltage	-			AC/DC 3.3 ± 0.05 V	
Active Voltage	-			3.0 - 4.0 V	
Consumption Current	-			0.38 A Typical (3.3 V)	

Pana-Max2

Model	Pana-Max2 M4
Hose Connection Type	ISO 9168 Type 2
Head Type	Standard Head
Rotation Speed	350,000 - 450,000 min ⁻¹
Bur Type	ISO 1797-1 Type 3 Ø1.59 - 1.60 mm Standard Bur
Chuck Type	Push Button
Chucking Length	10.2 mm
Max. Bur Length	21 mm
Max. Working Part Diameter	Ø2 mm
Recommended Drive Air Pressure	0.25 MPa (2.5 kgf/cm ²)
Operating Drive Air Pressure Range	0.20 - 0.25 MPa (2.0 - 2.5 kgf/cm ²)
Max. Air Consumption	< 66 NL/min (0.3±0.1 MPa)
Water Pressure	0.05 - 0.20 MPa (0.5 - 2.0 kgf/cm ²)

Model	Pana-Max2 M4
Max. Water Consumption	≥ 50 mL/min (0.20 MPa)
Chip Air Pressure	0.15 - 0.25 MPa (1.5 - 2.5 kgf/cm ²)
Rated Torque	≥ 0.05 cN·m (Recommended Max Drive Air Pressure)
Optic	-
Voltage	-
Active Voltage	-
Consumption Current	-

Pana-Max PLUS2 series

Model	PAP2-MU M4	PAP2-SU M4
Hose Connection Type	ISO 9168 Type 2	
Head Type	Miniature Head	Standard Head
Rotation Speed	380,000 - 450,000 min ⁻¹	325,000 - 410,000 min ⁻¹
Bur Type	ISO 1797-1 Type 3 Ø1.59 - 1.60 mm Short Shank Bur / Standard Bur	ISO 1797-1 Type 3 Ø1.59 - 1.60 mm Standard Bur
Chuck Type	Push Button	
Chucking Length	9 mm	10.5 mm
Max. Bur Length	21 mm	25 mm
Max. Working Part Diameter	Ø2mm	
Recommended Drive Air Pressure	0.25 MPa (2.5 kgf/cm ²)	
Operating Drive Air Pressure Range	0.25 - 0.30 MPa (2.5 - 3.0 kgf/cm ²)	
Max. Air Consumption	< 66 NL/min (0.3 ± 0.1 MPa)	

Model	PAP2-MU M4	PAP2-SU M4
Water Pressure	0.08 - 0.20MPa (0.8 - 2.0 kgf/cm ²)	
Max. Water Consumption	≥50 mL/min (0.20 MPa)	
Chip Air Pressure	0.10 - 0.20 MPa (1.0 - 2.0 kgf/cm ²)	
Rated Torque	≥0.05 cN·m (Recommended Max Drive Air Pressure)	
Optic	-	
Voltage	-	
Active Voltage	-	
Consumption Current	-	
Use Environment	Temperature : 10 - 40 °C Humidity : 30 - 75 %	
Transportation and Store Environment	Temperature : -10 - 50 °C Humidity : 10 - 85 % Pressure : 500 - 1,060 hPa	

M800L / M800 / M900L / M900 / S-Max pico For NSK Coupling
M800KL / M900KL / M900K / S-Max pico KL For KaVo® MULTiflex®, MULTiflex® LUX Coupling KaVo® and MULTiflex® are registered trademarks of KaVo Dental GmbH & Co.,Germany.
S-Max pico SL For Sirona® Quick Coupling Sirona® is a registered trademark of Sirona Dental Systems GmbH, Germany.
S-Max pico BLED For Bien Air® Unifix®, Unifix® L Coupling Bien Air® and Unifix® are registered trademarks of Bien-Air Dental S.A.,Switzerland.
M800WL / M900WL / M900W / S-Max pico WLED For W&H® Roto Quick® Coupling W&H® and Roto Quick® are registered trademarks of W&H Dentalwerk Bürmoos GmbH, Austria.

12-2 Symbol

Move in the direction indicated.

Move fully to the stop, in the direction indicated.

Prohibited / Process NOT permitted / Wrong

Do the step / Process permitted / Correct

State of being fixed.

Wear protective gloves

Wear a face mask

Wear eye protection

This product can be used with a washer-disinfector for thermal disinfection.

This product can be sterilized in a steam sterilizer at 135 °C.

Hand operated

Caution, hot surface

Refer to the operation manual

Caution, Consult operation instructions.

Date of manufacture

Conforms to CE European Directive of "Medical device directive 93/42/EEC."

Unique Device Identifier

Manufacturer.

GS1 DataMatrix for Unique Device Identifier.

Caution:U.S.Federal law restricts this device to sale by or on the order of a licensed physician.

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NSK America Corp. www.nskdental.com